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Warnings

1: The use of the device is limited only to medical staff specialized in dental implantology such as qualified surgeons or deontologists.

2: Once used, the device must be disposed in compliance with the regulations of the country in which it is used.

3: It is recommended to clean, disinfect and sterilize the devices before using them.

When tightening, use the recommended ratchet connection, handpiece driver and wrench. Ensure application of the recommended torques:

- Final abutments: tightening torque 35 Ncm

- Temporary abutments: tightening torque between 15-35 Ncm (do not tighten to 35 Ncm unless the implant is stable and fully osseointegrated).

Description

Prosthesis elements are single-use and made of medical Titanium. Before being packed, the product is exposed to a superficial decontamination treatment, wrapped inside plastic film bags and put in card board boxes. The labels, that are placed on both the bag and on the box, show code, description and technical characteristics of the product.

There are no particular rules for conservation: these MDs are intended to be used by non-profane users - professional studies of dental care specialists - which are suitable places to ensure proper conditions of conservation.

It is recommended that these MDs be manipulated in clean places to ensure hygienic conditions that are adequate to the type of activity carried out. For transport from the package to the patient's oral cavity use the compatible instrument and check the good status of the instrument itself.

Instructions for use

Prosthesis elements have been designed to be positioned in the implants to substitute dental crowns or to fix movable prosthesis.

Prosthesis elements are sold as non-sterile and they must be sterilized according to a validated procedure that must always be monitored.

Sterilization can be performed as follow:

-dry sterilization:

in oven at 160°C for 2 hrs;

-steam sterilization (autoclave) for USA:

seal the device in a bag and sterilize it by steam at 132°C (270°F) for 3 minutes;

make sure to comply the drying times of the device.

-steam sterilization (autoclave) outside the USA:

seal the device in a bag and sterilize it by steam at 132°C (270°F) for 3 minutes;

make sure to comply the drying times of the device.

Alternative for UK:

seal the device in a bag and sterilize it by steam at 134-135°C (273-275°F) for 3 minutes;

make sure to comply the drying times of the device.

Whereas steam sterilization stricter requirements are foreseen, local or national rules relative to the matter must be followed.

IMPIANTI DENTALI GIESSE
GIESSE DENTAL IMPLANTS
ELEMENTI PROTESICI/ 
PROSTHESIS ELEMENTS

16016 COGOLETO (GE) VIA MALUEA, 3 - PAD. 18
 TEL +39 010.9189127
 sales@giessetechnology.it www.giessetechnology.it

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The deontologist has the responsibility to verify its own sterilization program.
 For a good functioning, the accessories must be used exclusively for the preparation of the implant seat where the GIESSE implants will be implanted.
 The use of non-sterilized devices might induce to infection.

Precautions

The product is designed and built exclusively for the expected foreseen conditions; any other different use of it is to be considered inappropriate and dangerous.
 For good functionality, prosthesis elements must be used exclusively with the related GIESSE dental implant.

Surgical procedure:

First phase: take the healing screw off.

If the implant is submerged it will be necessary to incise the edge of the gum.

Second phase: insert the prosthesis element in the specific hole of the implant and fix it.

Third phase: do a dental impression of the prosthesis.

For more detailed information please see surgical protocol.

Warranties and limitations

GIESSE TECHNOLOGY Srl guarantees to the first buyer of this GIESSE prosthesis for dental implants that it was produced, packed and tested with reasonable care and that it has no manufacturing defects and not even in the materials for a period of two years starting from the date of the first purchase. GIESSE TECHNOLOGY Srl guarantees the products on condition that the packaging is not damaged or handled improperly. Defects consequent to the improper use of this product are not covered by this warranty. Moreover, the warranty is not applicable in case the accessory is used improperly or not in compliance with what is written in the instructions for use, impacts, falls, exposition to humidity or other environmental conditions or inadequate agents, or in any other case where the defect of the product is caused by a fact or act that might not be attributed, directly or indirectly, to GIESSE TECHNOLOGY Srl, as well as in case of removal, abrasion, alteration or the impossibility of identification of the product's lot number. Under no circumstances GIESSE TECHNOLOGY Srl will be believed to be responsible to the buyer for any direct or indirect damage related to the handling or use of the product, or for patients' choice. The responsibility in accordance with this warranty is limited to the replacement of the defective product or, exclusively at discretion of GIESSE TECHNOLOGY Srl, to the refund of the paid net amount.

	LOT CODE		WARNING, SEE THE INSTRUCTIONS FOR USE.
	CATALOGUE NUMBER (PRODUCT CODE)		WARNING (specific warnings or precautions)
	DO NOT REUSE (SINGLE USE)		LATEX FREE
	Do not use if the package is not intact		THESE DEVICES WITH THE CE MARK FULFILLS THE REQUIREMENTS OF THE MEDICAL DEVICES DIRECTIVE 93/42/EEC